

CONTAINER WITH FILTER large h 100 mm - red



Product Code	37110
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2811854
CND	S9099
EAN/UPC	8023279371109

Description

Standard container with two filters and non perforated bottom.

- Large model: 580 x 280 mm
- Height: 100 mm
- Non perforated bottom
- 2 filters
- Colour: red

Standard containers (supplied with 1 paper filter and filter holder).
All containers are autoclavable at maximum 134°C and 2.2 bar.

NOT FOR SALE IN GERMANY.