



|               |                |
|---------------|----------------|
| Product Code  | 37150          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|            |               |
|------------|---------------|
| RDM (NSIS) | 2811854       |
| CND        | S9099         |
| EAN/UPC    | 8023279371505 |

## Description

Barrier System container with valve and non perforated bottom.

This system is preferred when is requested continuous use of containers, it avoids filter replacement. Valve opens during sterilization to allow exit of steam.

Lock nail on bottom side made of chrome.

- Small model: 285 x 280 mm
- Height: 100 mm
- Non perforated bottom
- Barrier System with valve
- Colour: Grey

All containers are autoclavable at maximum 134°C and 2.2 bar.

NOT FOR SALE IN GERMANY.