

B.S. CONTAINER WITH VALVE medium h 135 mm - grey



Product Code	37201
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2811854
CND	S9099
EAN/UPC	8023279372014

Description

Barrier System container with valve and non perforated bottom.
This system is preferred when is requested continuous use of containers, it avoids filter replacement. Valve opens during sterilization to allow exit of steam.
Lock nail on bottom side made of chrome.

- Medium model: 465 x 280 mm
- Height: 135 mm
- Non perforated bottom
- Barrier System with valve
- Colour: Grey

All containers are autoclavable at maximum 134°C and 2.2 bar.

NOT FOR SALE IN GERMANY.