



Product Code	37212
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

EAN/UPC 8023279372120

## Description

Barrier System container with valve and non perforated bottom.

This system is preferred when is requested continuous use of containers, it avoids filter replacement. Valve opens during sterilization to allow exit of steam.

Lock nail on bottom side made of chrome.

- Medium model: 465 x 280 mm
- Height: 150 mm
- Non perforated bottom
- Barrier System with valve
- Colour: red

All containers are autoclavable at maximum 134°C and 2.2 bar.

NOT FOR SALE IN GERMANY.