

B.S. CONTAINER WITH VALVE large h 200 mm - grey



Product Code	37253
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2811854
CND	S9099
EAN/UPC	8023279372533
GMDN	13818

Description

Barrier System container with valve and non perforated bottom.
This system is preferred when is requested continuous use of containers, it avoids filter replacement. Valve opens during sterilization to allow exit of steam.
Lock nail on bottom side made of chrome.

- Large model: 580 x 280 mm
- Height: 200 mm
- Non perforated bottom
- Barrier System with valve
- Colour: grey

All containers are autoclavable at maximum 134°C and 2.2 bar.

NOT FOR SALE IN GERMANY.