

# B.S. CONTAINER WITH VALVE large h 100 mm - red



Product Code	37260
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

EAN/UPC 8023279372601

## Description

Barrier System container with valve and non perforated bottom.

This system is preferred when is requested continuous use of containers, it avoids filter replacement. Valve opens during sterilization to allow exit of steam.

Lock nail on bottom side made of chrome.

- Large model: 580 x 280 mm
- Height: 100 mm
- Non perforated bottom
- Barrier System with valve
- Colour: red

All containers are autoclavable at maximum 134°C and 2.2 bar.

NOT FOR SALE IN GERMANY.