

# B.S. CONTAINER WITH VALVE large h 100 mm - blue - perforated



Product Code	37280
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2811854
CND	S9099
EAN/UPC	8023279372809

## Description

Barrier System container with valve and perforated bottom.

This system is preferred when is requested continuous use of containers, it avoids filter replacement. Valve opens during sterilization to allow exit of steam.

Lock nail on bottom side made of chrome.

Perforated bottom allows much quicker sterilization process (half the time compared to not perforated).

- Large model: 580 x 280 mm
- Height: 100 mm
- Perforated bottom
- Barrier System with valve
- Colour: blue

All containers are autoclavable at maximum 134°C and 2.2 bar.

NOT FOR SALE IN GERMANY.