

# B.S. CONTAINER WITH VALVE large h 200 mm - blue - perforated



Product Code	37283
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

EAN/UPC 8023279372830

## Description

Barrier System container with valve and perforated bottom.

This system is preferred when is requested continuous use of containers, it avoids filter replacement. Valve opens during sterilization to allow exit of steam.

Lock nail on bottom side made of chrome.

Perforated bottom allows much quicker sterilization process (half the time compared to not perforated).

- Large model: 580 x 280 mm
- Height: 200 mm
- Perforated bottom
- Barrier System with valve
- Colour: blue

All containers are autoclavable at maximum 134°C and 2.2 bar.

NOT FOR SALE IN GERMANY.