

ETHICON VICRYL RAPID ABSORBABLE SUTURES - gauge 3/0 needle 17 mm - braided



Product Code	22373
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	III
UK-REP	No
CH-REP	No
RDM (NSIS)	2436959
CND	H0101010299
GTIN	30705031106073
GMDN	17471

Description

VICRYL RAPID BRAIDED - absorbable sterile surgical sutures

Complete absorption in 42 days

- Ethicon code: V2150G
- Gauge: 3/0
- Metric: 2.0
- Colour: undyed
- Circle: 1/2
- Needle length: 17 mm
- Needle: RB-1 PLus
- Suture length: 70 cm
- Minimum order: box of 12 pieces

Vicryl Rapid Braided indicated only for: use in superficial soft tissue approximation of the skin and mucosa, where only short term (7-10 days) wound support is required.