


# PLASTIC SECURITY SEAL - blue



Product Code	37330
--------------	-------

Unit of sale	box of 100 pcs
--------------	-------------------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

CH-REP	Not required
--------	--------------

EAN/UPC	8023279373301
---------	---------------

## Description

PLASTIC SECURITY SEALS - 2.7 cm LONG PIN

Security seal - blue

Container seals for steam sterilization processes.  
Conform to ISO 11140-Type 1.