



|               |                   |
|---------------|-------------------|
| Product Code  | 37340             |
| Unit of sale  | box of 100 pcs    |
| Minimum order | 1                 |
| Type          | No medical device |
| Class         |                   |

|        |    |
|--------|----|
| UK-REP | No |
| CH-REP | No |

|         |               |
|---------|---------------|
| EAN/UPC | 8023279373400 |
|---------|---------------|

## Description

PLASTIC SECURITY SEALS - 2.7 cm LONG PIN

Security seal - red

Container seals for steam sterilization processes.  
Conform to ISO 11140-Type 1.