

ETHICON VICRYL RAPID ABSORBABLE SUTURES - gauge 3/0 needle 19 mm - braided



Product Code 22375

Unit of sale box of 36 pcs

Minimum order 1

Type Medical device

Class III

UK-REP No

CH-REP No

RDM (NSIS) 2437028

CND H0101010299

GTIN 30705031150687

GMDN 17471

Description

VICRYL RAPID BRAIDED - absorbable sterile surgical sutures

Complete absorption in 42 days

- Ethicon code: VR2298
- Gauge: 3/0
- Metric: 2.0
- Colour: undyed
- Circle: 3/8
- Needle length: 19 mm
- Needle: FS-2
- Suture length: 75 cm
- Minimum order: box of 36 pieces

Vicryl Rapid Braided indicated only for: use in superficial soft tissue approximation of the skin and mucosa, where only short term (7-10 days) wound support is required.