


# PRE-CUT TOURNIQUET 45x2.5 cm - blue



Product Code	37510
Unit of sale	box of 250 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
--------	-----

RDM (NSIS)	2340978
CND	V9003
EAN/UPC	8023279375107

## Description

PRE-CUT TOURNIQUET 45x2.5 cm - blue - box of 250

Latex free disposable tape tourniquet made of TPE with smooth surface. Thickness 0.636 mm

Multilingual box: GB, FR, IT, ES, PT, DE, PL, GR, SE, HU, Arabic.