


PRE-CUT TOURNIQUET 45x2.5 cm - violet



Product Code	37512
Unit of sale	box of 250 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2340978
CND	V9003
EAN/UPC	8023279375121
PARAF	950833917

Description

PRE-CUT TOURNIQUET 45x2.5 cm - violet - box of 250

Latex free disposable tape tourniquet made of TPE with smooth surface. Thickness 0.636 mm

Multilingual box: GB, FR, IT, ES, PT, DE, PL, GR, SE, HU, Arabic.