

# TAPE TOURNIQUET - disposable, latex free, 25 pre-cut



Product Code	37515
Unit of sale	box of 25 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	No

RDM (NSIS)	2340978
CND	V9003
EAN/UPC	8023279375152

## Description

TAPE TOURNIQUET - latex free - disposable - box of 25

Latex free disposable tape tourniquet made of TPE with dot surface.

Each box contains 25 pre-cut tourniquets 45x2.5 cm (18"x1") with 0.636 mm thickness.

Multilingual box: GB, FR, IT, ES, PT, DE, PL, GR, SE, HU, Arabic.