

CHILD PROBE 3-20 kg Nellcor comp.- disposable



Product Code	33756
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	II B

RDM (NSIS)	1727075
CND	C900301
EAN/UPC	6945040104339
UMDNS	37808
PARAF	987472747

Description

Disposable SpO2 sensor

- For child 3-20 kg

Compatibility:

- NELLCOR (not compatible with Nellcor Oximax)

All sensors are supplied in sterile pack with 70 cm cable.

Technical Specifications

Range of SpO2 measurement: 70%~100%
Accuracy: 70%~100% ±2%; below 70%: unspecified
Range of pulse measurement: 30~250 bpm
Accuracy: ±2 bpm of 2% (select larger)
Optical Sensor:
red light (wavelength is 650~670 nm, 6.65 mW)
infrared (wavelength is 880~910 nm, 6.75 mW)