


# TERUMO SURFLO WINGED INFUSION SETS 21G - 0.8x19 mm - SV\*21NL30 - green - sterile



Product Code	37592
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2458168
CND	A01010202
GTIN	36927675502910

## Description

SURFLO™ WINGED INFUSION SETS - sterile

TERUMO® code: SV\*21NL30

Needle: 21G - 0.8×19 mm

Colour: green

Minimum order: box of 50

Surflo™ Winged Infusion Sets (Scalp Vein or SV set) are intended to introduce fluids into the vascular system after venipuncture. To ensure stability and grip during venipuncture, they are equipped with color-coded, interlocking wings. Thanks to the slightly angled wings (<math><3^\circ</math> in relation to the cannula), less pressure is put on the blood vessel when puncturing. They allow smooth skin penetration.

Not made with DEHP and have no natural rubber latex components.

Tube length: 30 cm.

Class IIa, Medical Device.

Multilingual box: GB, FR, IT, ES, PT, DE, PL, NL, HU, GR, BG, RO, FI, NO, SE, DK, CZ, SK, SI, RS, EE, LV, LT, TR, RU.