

TERUMO SURFLO WINGED INFUSION SETS 25G - 0.5x19 mm - SV*25NL30 - orange - sterile



| | |
|---------------|----------------|
| Product Code | 37595 |
| Unit of sale | box of 50 pcs |
| Minimum order | 1 |
| Type | Medical device |
| Class | II A |

| | |
|------------|----------------|
| RDM (NSIS) | 2458166 |
| CND | A01010202 |
| GTIN | 36927675502941 |

Description

SURFLO™ WINGED INFUSION SETS - sterile

TERUMO® code: SV*25NL30

Needle: 25G - 0.5×19 mm

Colour: orange

Minimum order: box of 50

Surflo™ Winged Infusion Sets (Scalp Vein or SV set) are intended to introduce fluids into the vascular system after venipuncture. To ensure stability and grip during venipuncture, they are equipped with color-coded, interlocking wings. Thanks to the slightly angled wings (<math><3^\circ</math> in relation to the cannula), less pressure is put on the blood vessel when puncturing. They allow smooth skin penetration.

Not made with DEHP and have no natural rubber latex components.

Tube length: 30 cm.

Class IIa, Medical Device.

Multilingual box: GB, FR, IT, ES, PT, DE, PL, NL, HU, GR, BG, RO, FI, NO, SE, DK, CZ, SK, SI, RS, EE, LV, LT, TR, RU.