

TERUMO AGANI HYPODERMIC NEEDLES 18Gx1 1/2" - 1.2x38mm - AN*1838R1 - pink - sterile



Product Code	37600
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	II A

CH-REP	Yes
--------	-----

RDM (NSIS)	795408
CND	A01010101
GTIN	36923033401338

Description

AGANI™ HYPODERMIC NEEDLES

- Terumo code: AN*1838R1
- Needle: 18Gx1 1/2" - 1.2x38 mm
- Colour: pink
- Minimum order: box of 100

Sterile hypodermic needles for single use, suitable for a wide range of applications. Special geometry of the sharp lancet point facilitates the penetration of the skin and helps to minimize the trauma. Silicone-coated cannula provides good gliding properties for a smooth penetration and prevents any knocking effect.

Agani hypodermic needles provide you and your patients with the best possible levels of care and comfort. Hub made of polypropylene, needle made of stainless steel. ISO colour-coded hub allows to recognize the diameter of the needle. Class II a. MDR 2017/745. Shelf life at production: 5 years. Individually packed in blister.

Multilingual box: GB, FR, IT, ES, PT, DE, PL, NL, HU, GR, BG, RO, FI, NO, SE, DK, CZ, SK, SI, RS, EE, LV, LT, TR, RU.