

# HYPODERMIC NEEDLE 26G 0.45x12.7 mm - sterile



Product Code	23761
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2357580
CND	A0101010102
EAN/UPC	8032764156268
PARAF	915924373

## Description

### STANDARD HYPODERMIC NEEDLES 26G - STERILE

Standard hypodermic needles, made in stainless steel AISI 304 with triple bevelling and totally non-toxic and non-pyrogenic. Each element of the needle is treated for painless injection and ensuring minimal trauma during penetration.

- Size: 26G - 0.45x13 mm
- Reference colour: brown
  
- Box of 100 pieces
- Individually packaged in blister
- Box and instructions: GB, IT, ES, PT.