

5-LEAD PATIENT CABLE



| | |
|--------------|-------|
| Product Code | 33763 |
|--------------|-------|

| | |
|--------------|------|
| Unit of sale | 1 pc |
|--------------|------|

| | |
|---------------|---|
| Minimum order | 1 |
|---------------|---|

| | |
|------|--|
| Type | |
|------|--|

| | |
|-------|--|
| Class | |
|-------|--|

| | |
|---------|---------------|
| EAN/UPC | 8809276946300 |
|---------|---------------|

Description

ECG accessories for monitor BM3 and BM5 (codes: 33718,19,20,21,22,23,24,80,81).

- 5-lead cable for kit 33712 - spare
- New model (after 6/2006)