

# TERUMO SYRINGES W/O NEEDLE 10 ml - Concentric Luer Lock - SS\*10LE1 - sterile



Product Code	37634
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	II A

CH-REP	Yes
RDM (NSIS)	2431057
CND	A020102010202
GTIN	34987350702648

## Description

- TERUMO® SYRINGES WITHOUT NEEDLE - STERILE
- CONCENTRIC LUER LOCK
- TERUMO® code: SS\*10LE1 (Compatible with Terumo infusion pump)
- Capacity: 10 ml
- Minimum order: box of 100

Sterile 3-part syringes for single use, with or without pre-connected sterile hypodermic needle.  
Syringes are intended for the aspiration of fluid or blood, or for the injection of fluids immediately after filling into the human body.  
Hub is made of polypropylene and master batch, while cannula is made of stainless steel.  
Latex free.  
Transparent barrel provides easy recognition of air bubbles and the clear printing allows for precise dosage management.

Class IIa, Medical Device following MDR 2017/745.  
Shelf life at production: 5 years.  
Multilingual box: GB, FR, IT, ES, PT, DE, PL, NL, HU, GR, BG, RO, FI, NO, SE, DK, CZ, SK, SI, RS, EE, LV, LT, TR, RU.