

# TERUMO SYRINGES W/O NEEDLE 30 ml - Concentric Luer Lock - SS+30L1MP - sterile



Product Code	37637
Unit of sale	box of 25 pcs
Minimum order	1
Type	Medical device
Class	II A

CH-REP	Yes
GTIN	34806017513369

## Description

- TERUMO® SYRINGES WITHOUT NEEDLE - STERILE
- CONCENTRIC LUER LOCK
  - TERUMO® code: SS+30L1MP (Compatible with Terumo infusion pump)
  - Capacity: 30 ml
  - Minimum order: box of 25

Sterile 3-part syringes for single use, with or without pre-connected sterile hypodermic needle. Syringes are intended for the aspiration of fluid or blood, or for the injection of fluids immediately after filling into the human body. Hub is made of polypropylene and master batch, while cannula is made of stainless steel. Latex free. Transparent barrel provides easy recognition of air bubbles and the clear printing allows for precise dosage management.

Class IIa, Medical Device following MDR 2017/745.

Shelf life at production: 5 years.

Multilingual box: GB, FR, IT, ES, PT, DE, PL, NL, HU, GR, BG, RO, FI, NO, SE, DK, CZ, SK, SI, RS, EE, LV, LT, TR, RU.