

# TERUMO SYRINGES W/O NEEDLE 1 ml Tuberculin - Concentric Luer Slip - MDSS01SE - sterile



|               |                |
|---------------|----------------|
| Product Code  | 37640          |
| Unit of sale  | box of 100 pcs |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | II A           |

|        |     |
|--------|-----|
| CH-REP | Yes |
|--------|-----|

## Description

- TERUMO® SYRINGES WITHOUT NEEDLE - STERILE
- CONCENTRIC LUER SLIP
  - TERUMO® code: MDSS01SE
  - Capacity: 1 ml - Tuberculin
  - Minimum order: box of 100

Sterile 3-part syringes for single use, with or without pre-connected sterile hypodermic needle. Syringes are intended for the aspiration of fluid or blood, or for the injection of fluids immediately after filling into the human body. Hub is made of polypropylene and master batch, while cannula is made of stainless steel. Latex free. Transparent barrel provides easy recognition of air bubbles and the clear printing allows for precise dosage management.

Class IIa, Medical Device following MDR 2017/745.

Shelf life at production: 5 years.

Multilingual box: GB, FR, IT, ES, PT, DE, PL, NL, HU, GR, BG, RO, FI, NO, SE, DK, CZ, SK, SI, RS, EE, LV, LT , TR, RU.