

TERUMO SYRINGES W/O NEEDLE 20 ml - Eccentric Luer Slip - MDSS20ESE - sterile



Product Code	37646
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	II A

CH-REP	Yes
--------	-----

RDM (NSIS)	2059606
CND	A020102010202
GTIN	34806017500956

Description

- TERUMO® SYRINGES WITHOUT NEEDLE - STERILE
- ECCENTRIC LUER SLIP
- TERUMO® code: MDSS20ESE
- Capacity: 20 ml
- Minimum order: box of 50

Sterile 3-part syringes for single use, with or without pre-connected sterile hypodermic needle. Syringes are intended for the aspiration of fluid or blood, or for the injection of fluids immediately after filling into the human body. Hub is made of polypropylene and master batch, while cannula is made of stainless steel. Latex free. Transparent barrel provides easy recognition of air bubbles and the clear printing allows for precise dosage management.

Class IIa, Medical Device following MDR 2017/745.
Shelf life at production: 5 years.
Multilingual box: GB, FR, IT, ES, PT, DE, PL, NL, HU, GR, BG, RO, FI, NO, SE, DK, CZ, SK, SI, RS, EE, LV, LT, TR, RU.