

TERUMO SYRINGES W/O NEEDLE 50 ml - Catheter tip - SS+50C1 - sterile



Product Code	37648
Unit of sale	box of 25 pcs
Minimum order	1
Type	Medical device
Class	II A

CH-REP	Yes
--------	-----

GTIN	34806017507849
------	----------------

Description

TERUMO® SYRINGES WITHOUT NEEDLE - STERILE

- CATHETER TIP
- TERUMO® code: SS+50C1
- Capacity: 50 ml
- Minimum order: box of 25

Sterile 3-part syringes for single use, with or without pre-connected sterile hypodermic needle.

Syringes are intended for the aspiration of fluid or blood, or for the injection of fluids immediately after filling into the human body.

Hub is made of polypropylene and master batch, while cannula is made of stainless steel.

Latex free.

Transparent barrel provides easy recognition of air bubbles and the clear printing allows for precise dosage management.

Class IIa, Medical Device following MDR 2017/745.

Shelf life at production: 5 years.

Multilingual box: GB, FR, IT, ES, PT, DE, PL, NL, HU, GR, BG, RO, FI, NO, SE, DK, CZ, SK, SI, RS, EE, LV, LT, TR, RU.