

# TERUMO SYRINGES WITH NEEDLE 3 ml - 21Gx1/2" - MDSS03S2138E - sterile



Product Code	37654
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	Yes
RDM (NSIS)	2303174
CND	A020102010201
GTIN	36970114923784

## Description

- TERUMO® SYRINGES WITH NEEDLE - STERILE
- TERUMO® code: MDSS03S2138E
  - Capacity: 3 ml
  - Assembled needle: 21Gx1 1/2" - 0.8x38 mm
  - Minimum order: box of 100

Sterile 3-part syringes for single use, with or without pre-connected sterile hypodermic needle. Syringes are intended for the aspiration of fluid or blood, or for the injection of fluids immediately after filling into the human body. Hub is made of polypropylene and master batch, while cannula is made of stainless steel. Latex free. Transparent barrel provides easy recognition of air bubbles and the clear printing allows for precise dosage management.

Class IIa, Medical Device following MDR 2017/745.

Shelf life at production: 5 years.

Multilingual box: GB, FR, IT, ES, PT, DE, PL, NL, HU, GR, BG, RO, FI, NO, SE, DK, CZ, SK, SI, RS, EE, LV, LT, TR, RU.