

# TERUMO SYRINGES WITH NEEDLE 1 ml Tuberculin - 26Gx1/2" - MDSS01S2613E - sterile



Product Code	37662
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	II A

CH-REP	Yes
--------	-----

RDM (NSIS)	2059438
CND	A020102010201
PARAF	983742329
GTIN	36970114923043

## Description

- TERUMO® SYRINGES WITH NEEDLE - STERILE  
- TUBERCULIN  
- TERUMO® code: MDSS01S2613E  
- Capacity: 1 ml  
- Assembled needle: 26Gx1/2" - 0.45x13mm  
- Minimum order: box of 100

Sterile 3-part syringes for single use, with or without pre-connected sterile hypodermic needle. Syringes are intended for the aspiration of fluid or blood, or for the injection of fluids immediately after filling into the human body. Hub is made of polypropylene and master batch, while cannula is made of stainless steel.  
Latex free.

Transparent barrel provides easy recognition of air bubbles and the clear printing allows for precise dosage management.

Class IIa, Medical Device following MDR 2017/745.

Shelf life at production: 5 years.

Multilingual box: GB, FR, IT, ES, PT, DE, PL, NL, HU, GR, BG, RO, FI, NO, SE, DK, CZ, SK, SI, RS, EE, LV, LT, TR, RU.