

INSTRUMENT TRAY - MESH BASE 300x200x52 mm - plastic



Product Code	37691
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

GTIN 15060098131293

Description

PLASTIC HOLLOWARE - Autoclavable hundreds of times

Instrument tray - mesh base and sides
Size mm: 300x250x52
Graduation: -

For the protection and delivery of sterile instruments to operating rooms. For the administration of procedures, holding cleansing solutions and the administration of medicines and fluids.

4 good reasons to choose plastic line:

- 50% lighter than stainless steel
- autoclavable: 1000 cycles guaranteed
- attractive modern design
- cheaper than equivalent stainless steel

Made of polypropylene material, easy to clean, resistant to most healthcare chemicals. Machine, ultrasonic or manual cleaning are permissible. Autoclave sterilization at 134.6°C (+2°C) for a minimum of 2 minutes. EO Ethylene Oxide or Hydrogen Peroxide sterilization both at max. 65°C.

Made in the UK.