

GALLIPOT 240 ml - plastic



Product Code	37702
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2740271
CND	V0402
EAN/UPC	8023279377026
GTIN	08023279377026
GMDN	42893

Description

GALLIPOT 63x42 mm - 90 ml

Economic line of reusable plasticware products.
Made in polypropylene, easy to clean and resistant to most healthcare chemicals.

Autoclavable at 121°C (15 min), 126°C (10 min) or 134°C (3 min).