

# KIDNEY TRAY 6" 155x75 mm - plastic



Product Code	37705
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2739661
CND	V0402
EAN/UPC	8023279377057
GTIN	08023279377057
GMDN	42893

## Description

KIDNEY TRAY 6" 155 x 75 mm - 250 ml

Economic line of reusable plasticware products.  
Made of polypropylene, easy to clean and resistant to most healthcare chemicals.

Autoclavable at 121°C (15 min), 126°C (10 min) or 134°C (3 min).