

KIDNEY TRAY 12" 306x140 mm - plastic



Product Code	37708
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2739661
CND	V0402
EAN/UPC	8023279377088
GTIN	08023279377088
GMDN	42893

Description

KIDNEY TRAY 12" 306 x 140 mm - 1500 ml

Economic line of reusable plasticware products.
Made in polypropylene, easy to clean and resistant to most healthcare chemicals.

Autoclavable at 121°C (15 min), 126°C (10 min) or 134°C (3 min).