

# MEASURING JUG 500 ml - plastic



Product Code	37710
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2800782
CND	V0899
EAN/UPC	8023279377101
PARAF	987471354
GTIN	08023279377101
GMDN	66310

## Description

Measuring jug 500 ml (raised graduation for easy reading).

Economic line of reusable plasticware products.  
Made in polypropylene, easy to clean and resistant to most healthcare chemicals.

Autoclavable at 121°C (15 min), 126°C (10 min) or 134°C (3 min).