

# INSTRUMENT TRAY WITH COVER 220x150x70 mm - plastic



Product Code	37718
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2740094
CND	V0402
EAN/UPC	8023279377187
GTIN	08023279377187
GMDN	12143

## Description

INSTRUMENT TRAY 220x150xh 70 mm

Provided with cover which fits into groove of the tray to avoid minimum spillage of material.

Economic line of reusable plasticware products.  
Made in polypropylene, easy to clean and resistant to most healthcare chemicals.

Autoclavable at 121°C (15 min), 126°C (10 min) or 134°C (3 min).