

S/S GOLDEN KIDNEY DISH 6"



Product Code	37750
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	No

RDM (NSIS)	2804713
CND	Z129099
EAN/UPC	8023279377507

Description

S/S AISI 202 "GOLDEN" KIDNEY DISHES

Tray 6"- 150 mm
Capacity: 200 ml

Golden kidney dishes easy to clean an sterilize. Nickel and chromium free, makes it hypoallergenic. Coating prevents corrosion.