



Product Code	37752
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	No
RDM (NSIS)	2804713
CND	Z129099
EAN/UPC	8023279377521

## Description

S/S AISI 202 "GOLDEN" KIDNEY DISHES

Tray 10"- 250 mm  
Capacity: 800 ml

Golden kidney dishes easy to clean and sterilize. Nickel and chromium free, makes it hypoallergenic. Coating prevents corrosion.