

# S/S GOLDEN KIDNEY DISH 12"



Product Code	37753
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	No

RDM (NSIS)	2804713
CND	Z129099
EAN/UPC	8023279377538

## Description

S/S AISI 202 "GOLDEN" KIDNEY DISHES

Tray 12" - 300 mm  
Capacity: 1500 ml

Golden kidney dishes easy to clean and sterilize. Nickel and chromium free, makes it hypoallergenic. Coating prevents corrosion.