



Product Code	37815
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

EAN/UPC	8023279378153
---------	---------------

## Description

### S/S KICK BUCKET - 12l

Kick bucket made of durable and rust-resistant stainless steel. Equipped with ergonomic anti-slip handle and 4 rubber swivelling castors, Ø 50 mm.

Easy to clean.

Overall size: Ø 350xh 340 mm.

Bucket size: Ø 300xh 260 mm.