

HYPODERMIC SAFETY NEEDLE 18G 1.2x38 mm - sterile



Product Code	23790
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2352379
CND	A0101010101
EAN/UPC	8034055916421

Description

HYPODERMIC SAFETY NEEDLES - STERILE

A line of hypodermic needles and scalp vein sets ensuring minimal trauma during penetration thanks to the triple bevelling.
Each element of the needle is treated for painless injection.
Needles are in stainless steel AISI 304 and totally non-toxic and non-pyrogenic.

- Size: 18G - 1.2x38 mm
- Reference colour: pink
- Box of 100 pieces
- Individually packaged in blister
- Box and instructions: GB, IT, ES, PT.