

TEST TUBE 12x75 mm - 5 ml - cylindrical, no rim



Product Code	37902
Unit of sale	box of 4000 pcs
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279379020
PARAF	977816002

Description

Polystyrene Test Tubes - graduated

- 5 ml cylindrical, no rim
- Dimensions: 12 x 75 mm

Box of 4,000 pieces (16 bag of 250 pcs.)

In order to fulfill every laboratory operator's necessity, a wide range of disposable test tubes, microtubes and pipettes is available in different sizes and materials. High quality raw materials, transparency and resistance are the main technical features that characterize these products.