

TEST TUBE 16x100 mm - 10 ml - cylindrical, with rim



Product Code	37912
Unit of sale	box of 2000 pcs
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279379129
PARAF	977816026

Description

Polystyrene Test Tubes - graduated

- 10 ml cylindrical, with rim
- Dimensions: 16x100 mm

Box of 2,000 pieces

In order to fulfill every laboratory operator's necessity, a wide range of disposable test tubes, microtubes and pipettes is available in different sizes and materials. High quality raw materials, transparency and resistance are the main technical features that characterize these products.