

TEST TUBE 5 ml - cylindrical - 12x86 mm - polypropylene - sterile



Product Code	37918
Unit of sale	box of 3000 pcs
Minimum order	1
Type	IVD - Professional Use
Class	

UK-REP	No
CH-REP	No

RDM (NSIS)	1911577
CND	W050301020102
EAN/UPC	8023279379181

Description

Test tube 5 ml - cylindrical - 12x86 mm

Sterile test tubes
Polypropylene test tubes

Versatile test tubes with red caps, which indicate that tubes are sterile.
All tubes are provided with writing area.
Beta rays sterilization.