

# TEST TUBE 10 ml - cylindrical - 16x100 mm - polypropylene - sterile



Product Code	37919
Unit of sale	box of 1600 pcs
Minimum order	1
Type	IVD - Professional Use
Class	

UK-REP	No
CH-REP	No

RDM (NSIS)	1498768
CND	W050301020102
EAN/UPC	8023279379198

## Description

Test tube 10 ml - cylindrical - 16x100 mm

Sterile test tubes  
Polypropylene test tubes

Versatile test tubes with red caps, which indicate that tubes are sterile.  
All tubes are provided with writing area.  
Beta rays sterilization.