

TEST TUBE 5 ml - cylindrical - 12x86 mm - polypropylene
- sterile



Product Code	37920
Unit of sale	box of 1500 pcs
Minimum order	1
Type	IVD - Professional Use
Class	

UK-REP	No
CH-REP	No

RDM (NSIS)	1911577
CND	W050301020102
EAN/UPC	8023279379204

Description

Polypropylene test tubes - sterile

- Test tube 5 ml
- cylindrical
- 12x86 mm
- box of 1,500

Versatile test tubes with red caps, which indicate that tubes are sterile.
All tubes are provided with writing area.
Beta rays sterilization.