

STOPPER for tube diam.16 mm



Product Code	37932
Unit of sale	box of 1000 pcs
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279379327
PARAF	977816065

Description

Stopper for tubes Ø 16 mm

Stoppers in polyethylene, neutral colour.

In order to fulfill every laboratory operator's necessity, a wide range of disposable test tubes, microtubes and pipettes is available in different sizes and materials. High quality raw materials, transparency and resistance are the main technical features that characterize these products.