

MICROTUBE 1.5 ml - conical



Product Code	37941
Unit of sale	box of 10000 pcs
Minimum order	1
Type	IVD - Professional Use
Class	A

UK-REP	No
CH-REP	No

RDM (NSIS)	2366874
CND	W050301020102
EAN/UPC	8023279379419

Description

Microtube 1.5 ml
box of 10,000 (20 pouches of 500)

Graduated microtubes made of polypropylene.
Flat cap, one handed usage, neutral colour, frosted area and perfect closure.