

MICROTUBE 2.0 ml - conical



Product Code	37942
Unit of sale	box of 8000 pcs
Minimum order	1
Type	IVD - Professional Use
Class	A

UK-REP	No
CH-REP	No

RDM (NSIS)	2373731
CND	W050301020102
EAN/UPC	8023279379426

Description

Microtube 2 ml
box of 8,000 (16 pouches of 500)

Graduated microtubes made of polypropylene.
Flat cap, one handed usage, neutral colour, frosted area and perfect closure.