

CENTRIFUGUE TEST TUBE 15 ml - conical - 17x120mm - non sterile



Product Code	37960
Unit of sale	box of 1700 pcs
Minimum order	1
Type	IVD - Professional Use
Class	A

UK-REP	No
CH-REP	No

RDM (NSIS)	2372087
CND	W050301020102
EAN/UPC	8023279379600

Description

Centrifuges test tubes in polypropylene - graduated with cap

Test tube 15 ml - conical - 17x120 mm
Box of 1,700
Non sterile

All test tubes have a writing area and screw cap.