

# CENTRIFUGUE TEST TUBE 50 ml - conical - 30x115mm - non sterile



Product Code	37961
Unit of sale	box of 500 pcs
Minimum order	1
Type	IVD - Professional Use
Class	A

UK-REP	No
CH-REP	No

RDM (NSIS)	2372093
CND	W050301020102
EAN/UPC	8023279379617

## Description

Centrifuges test tubes in polypropylene - graduated with cap

Test tube 50 ml - conical - 30x115 mm  
Box of 500  
Non sterile

All test tubes have a writing area and screw cap.