

CENTRIFUGUE TEST TUBE 50 ml - conical with base -
30x115 mm - non sterile



Product Code	37962
Unit of sale	box of 500 pcs
Minimum order	1
Type	IVD - Professional Use
Class	A

UK-REP	No
CH-REP	No

RDM (NSIS)	2372089
CND	W050301020102
EAN/UPC	8023279379624

Description

Centrifuges test tubes in polypropylene - graduated with cap

Test tube 50 ml - conical with base (size: 30x115 mm)
Box of 500
Non sterile

All test tubes have a writing area and screw cap.